

**LABORATORY FTIR / DATA SHEET**
**SF700Z**

# Fourier Transform Infrared Spectrometer

A high-energy, dynamically aligned FTIR platform engineered for dependable identification, advanced sampling and research workflows.

**Designed for confident routine analysis**

SF700Z combines a 24 W stabilized infrared source, four-quadrant dynamic collimation, high-sensitivity DLATGS detection and direct temperature/humidity monitoring.

- DSP-controlled non-orthogonal Michelson interferometer
- Four-quadrant dynamic collimation technology
- Digital temperature/humidity display with alarm
- ATR, reflection and infrared-microscope expansion



**520 × 430 × 250 mm**  
DIMENSIONS

**19 kg**  
WEIGHT

**< 30 W**  
NORMAL POWER

**7800–350**  
SPECTRAL RANGE,  
CM<sup>-1</sup>

**Key Advantages**

- DSP-controlled dynamic alignment supports optical stability and energy optimization.
- 24 W stabilized SiC source provides strong energy for demanding workflows.
- Sealed cavity and reusable 304 stainless-steel desiccant manage moisture.
- Flexible chamber supports transmission, ATR, reflection and microscopy accessories.

**Application Scenarios**

- Biopharma, materials science and petrochemical analysis.
- Food safety, environmental and polymer workflows.
- Coatings, surfaces, teaching and research.
- Identification across solids, powders, liquids and surfaces.

**MEASUREMENT PRINCIPLE**

## FTIR Platform

SF700Z combines DSP-controlled interferometry, four-quadrant dynamic collimation, a fully sealed optical cavity and guided software for dependable routine, industrial and research analysis.

|                              |                                    |                                 |                              |                                |
|------------------------------|------------------------------------|---------------------------------|------------------------------|--------------------------------|
| <b>01</b><br>24 W SiC source | <b>02</b><br>DSP<br>interferometer | <b>03</b><br>Sample / accessory | <b>04</b><br>DLATGS detector | <b>05</b><br>Fourier transform |
|------------------------------|------------------------------------|---------------------------------|------------------------------|--------------------------------|

**Engineered for stability, energy efficiency and easy maintenance**
**01 DSP-controlled interferometer**

A modified non-orthogonal Michelson interferometer uses DSP-controlled electromagnetic drive, continuous dynamic adjustment and automatic energy optimization.

**02 Dynamic four-quadrant collimation**

Laser-based four-quadrant alignment compensates moving-mirror mechanical deviation and optical errors caused by gravity or temperature change.

|                                                                                                                                                                                             |                                                                                                                                                                                                  |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>03 High-energy infrared source</b><br>A 24 W electronically stabilized silicon-carbide source and internal reflector provide strong energy for ATR and infrared-microscopy applications. | <b>04 Sealed cavity &amp; desiccant</b><br>The sealed optical platform and reusable 304 stainless-steel boxed desiccant manage moisture; the desiccant is accessible without removing the cover. |
| <b>05 Digital environment monitor</b><br>An industrial-grade temperature/humidity module provides direct digital readout and a humidity alarm for clear cavity-status awareness.            | <b>06 Flexible expansion</b><br>The large sample chamber supports transmission, ATR, diffuse, planar and external reflection, plus infrared-microscope expansion.                                |

## Bench-ready system at a glance

|                                 |                                           |                                     |                                   |
|---------------------------------|-------------------------------------------|-------------------------------------|-----------------------------------|
| <b>24-bit</b><br>A/D CONVERSION | <b>500 kHz</b><br>SAMPLING<br>ELECTRONICS | <b>24 W</b><br>STABILIZED IR SOURCE | <b>&lt; 30 W</b><br>TYPICAL POWER |
|---------------------------------|-------------------------------------------|-------------------------------------|-----------------------------------|

## TECHNICAL SPECIFICATIONS

# Performance Data

The following specifications are compiled from the supplied SF700Z technical document. Configuration-specific values should be confirmed at quotation.

## Optical performance

| Item                  | Specification                                                                                              |
|-----------------------|------------------------------------------------------------------------------------------------------------|
| Spectral range        | 7800–350 cm <sup>-1</sup> (approx. 1.28–28.57 μm)                                                          |
| Spectral resolution   | Better than 0.85 cm <sup>-1</sup> ; 0.5 cm <sup>-1</sup> optional                                          |
| Wavenumber precision  | ≤ 0.01 cm <sup>-1</sup>                                                                                    |
| Signal-to-noise ratio | Better than 45,000:1 (peak-to-peak, 4 cm <sup>-1</sup> , 1-minute acquisition)                             |
| Interferometer        | Modified non-orthogonal Michelson; DSP-controlled electromagnetic drive with continuous dynamic adjustment |
| Optical alignment     | Four-quadrant dynamic collimation with high-reflectivity gold-coated corner-reflector optics               |
| Beam splitter         | Imported germanium-coated KBr substrate with anti-fog moisture protection                                  |
| Detector              | Imported high-sensitivity DLATGS detector with moisture barrier                                            |
| Infrared source       | Imported 24 W electronically stabilized, high-power silicon-carbide rod source                             |

## System and electrical

| Item                | Specification                                                                     |
|---------------------|-----------------------------------------------------------------------------------|
| Laser system        | Imported low-power VCSEL, bandwidth better than 1 nm; high-stability He-Ne option |
| Digital electronics | SuperTect full-digital architecture; 24-bit, up to 500 kHz A/D conversion         |
| Environment monitor | Industrial-grade digital temperature/humidity display with humidity alarm         |

| Item                | Specification                                                                                     |
|---------------------|---------------------------------------------------------------------------------------------------|
| Moisture management | Fully sealed cavity; reusable 304 stainless-steel boxed desiccant with quick external access      |
| Data interface      | USB 2.0; compatible with USB 3.0 host ports                                                       |
| Sample chamber      | Large sliding-door compartment; supports transmission, ATR, reflection and microscopy accessories |
| Power consumption   | Less than 30 W during normal operation                                                            |
| Input power         | 100–240 VAC, 50/60 Hz                                                                             |
| Dimensions / weight | 520 × 430 × 250 mm / 19 kg                                                                        |
| Software platform   | Windows PC; supported operating-system versions depend on software release                        |

## Recommended environment

| Item              | Specification                                                                                        |
|-------------------|------------------------------------------------------------------------------------------------------|
| Laboratory        | Dry, clean, well ventilated and well sealed; protect the instrument from direct sunlight and dust    |
| Humidity control  | Use the built-in temperature/humidity display and humidity alarm; provide dehumidification as needed |
| Electrical        | Stable, grounded 100–240 VAC, 50/60 Hz supply                                                        |
| Sample area       | Keep the sample chamber, optical windows and accessories clean and dry                               |
| Long-term storage | Keep the optical cavity sealed and maintain the reusable desiccant in serviceable condition          |

## CONFIGURATION NOTE

Beam splitter, spectral libraries, sampling accessories and current software/OS compatibility are configuration dependent. Confirm the ordered configuration before purchase.

## Configuration options

| Option               | Source-supported configuration                            |
|----------------------|-----------------------------------------------------------|
| Resolution           | 0.5 cm <sup>-1</sup> high-resolution configuration        |
| Protected optics     | KRS-5 and ZnSe moisture-resistant material configurations |
| Sampling accessories | ATR plus diffuse, planar and external reflection          |
| Infrared microscopy  | Microscope expansion for localized analysis               |

## SOFTWARE AND DATA

# Intelligent Workflow

## From acquisition to confident identification

The infrared workstation combines guided acquisition, full-process preview, spectrum processing, multi-mode library search and reporting in one approachable workflow.

## Live acquisition preview

Preview spectra in real time while monitoring instrument status and the built-in temperature/humidity display.

|                                                                                                                                                                           |                                                                                                                                                                                      |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>1 Acquire and monitor</b><br>Guided acquisition and full-process preview monitor both background and sample spectra while instrument and cavity status remain visible. | <b>2 Process and interpret</b><br>Peak labeling, spectrum/constant arithmetic, spectrum/spectrum arithmetic, ATR correction, smoothing and derivative tools support interpretation.  |
| <b>3 Search and compare</b><br>Search by substance name, measured spectrum or file; compare files and identify samples against supplied or user-maintained libraries.     | <b>4 Calibrate and report</b><br>One-click intelligent laser-frequency correction supports standardized workflows, while result reporting and archiving keep routine work organized. |

## Spectral libraries, calibration and user knowledge

| Capability                 | Source-supported content                                                               |
|----------------------------|----------------------------------------------------------------------------------------|
| Standard reference spectra | More than 1,700 common substances in the built-in library                              |
| Optional libraries         | National pharmacopoeia and pharmaceutical packaging-material libraries                 |
| User-defined libraries     | Add, organize and maintain customer spectra as the knowledge base grows                |
| Laser correction           | One-click correction supports workflows aligned with GB/T 21186-2007 and JJF 1319-2011 |

## SAMPLING AND APPLICATIONS

# Flexible Configuration


## Three practical routes from sample to spectrum

Configure SF700Z for classic transmission, rapid ATR, reflection or infrared-microscopy workflows across solids, powders, liquids, films and surfaces.



### KBr pellet transmission POWDERS AND SOLIDS

Grind approximately 1 mg of dry sample with approximately 200 mg of dry spectroscopic-grade KBr, press at about 15 MPa for 30–60 seconds, then measure the transparent or translucent pellet.



### Demountable liquid cell FLOWING AND LOW-VISCOSITY LIQUIDS

Place a drop between compatible infrared windows with the required spacer, assemble gradually and measure. Clean the windows with a suitable solvent immediately after use and protect them from moisture.



### ATR measurement SOLIDS, POWDERS, LIQUIDS AND SURFACES

Place the sample in firm contact with a ZnSe or diamond ATR crystal and acquire the spectrum directly. ATR simplifies difficult-to-grind, insoluble and non-destructive surface measurements.

## Representative application areas

|                                                                                            |                                                                                                   |
|--------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------|
| <b>Biopharma</b><br>Identity and comparison of active materials and packaging              | <b>Materials science</b><br>Structure and composition studies across engineered materials         |
| <b>Petrochemicals</b><br>Qualitative screening of fuels, oils, additives and intermediates | <b>Food safety</b><br>Ingredient, additive and authenticity-oriented analysis                     |
| <b>Environmental</b><br>Screening of collected materials and contamination samples         | <b>Polymers &amp; rubber</b><br>Direct ATR analysis of difficult-to-grind or insoluble samples    |
| <b>Coatings &amp; surfaces</b><br>ATR, planar, diffuse and external-reflection workflows   | <b>Teaching &amp; research</b><br>Molecular identification, interpretation and method development |

## Common accessory selection

| Route               | Typical accessory set                                                           |
|---------------------|---------------------------------------------------------------------------------|
| Pellet preparation  | Press, pellet die, agate mortar, spectroscopic-grade KBr and drying oven        |
| Liquid transmission | Demountable liquid cell, infrared windows and thickness spacers                 |
| Advanced modes      | ZnSe or diamond ATR, diffuse/planar/external reflection and infrared microscope |

## VERIFICATION AND INSTALLATION

# Documented Performance

## Performance verification and lab guidance

The supplied source document includes peak-position checks and operating recommendations. The values below are reproduced as representative verification data, not as an individual instrument certificate.

## Peak-position verification

| Acceptance band              | Reference peak         | Measured peak             | Deviation               | Result |
|------------------------------|------------------------|---------------------------|-------------------------|--------|
| $\leq \pm 5 \text{ cm}^{-1}$ | $3027 \text{ cm}^{-1}$ | $3026.49 \text{ cm}^{-1}$ | $-0.51 \text{ cm}^{-1}$ | Pass   |
| $\leq \pm 5 \text{ cm}^{-1}$ | $2851 \text{ cm}^{-1}$ | $2850.44 \text{ cm}^{-1}$ | $-0.56 \text{ cm}^{-1}$ | Pass   |
| $\leq \pm 1 \text{ cm}^{-1}$ | $1601 \text{ cm}^{-1}$ | $1601.17 \text{ cm}^{-1}$ | $+0.17 \text{ cm}^{-1}$ | Pass   |
| $\leq \pm 1 \text{ cm}^{-1}$ | $1028 \text{ cm}^{-1}$ | $1028.63 \text{ cm}^{-1}$ | $+0.63 \text{ cm}^{-1}$ | Pass   |
| $\leq \pm 1 \text{ cm}^{-1}$ | $907 \text{ cm}^{-1}$  | $906.93 \text{ cm}^{-1}$  | $-0.07 \text{ cm}^{-1}$ | Pass   |

## RESOLUTION CHECK

The supplied verification document reports that seven peaks in the  $3110\text{--}2850 \text{ cm}^{-1}$  region were clearly resolved. SF700Z is specified at better than  $0.85 \text{ cm}^{-1}$  spectral resolution, with a  $0.5 \text{ cm}^{-1}$  option.

## Installation and care

|                                                                                                                                                                                                                                                                                                                                                                                            |                                                                                                                                                                                                                                                                                                                                                                                                                                      |                                                                                                                                                                                                                                                                                                                                                                                                              |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Laboratory setup</b></p> <ul style="list-style-type: none"> <li>• Use a dry, ventilated, clean room with good door and window sealing.</li> <li>• Protect the instrument from direct sunlight, dust and corrosive vapors.</li> <li>• Provide air conditioning and dehumidification according to local conditions.</li> <li>• Use a grounded 100–240 VAC, 50/60 Hz supply.</li> </ul> | <p><b>Moisture management</b></p> <ul style="list-style-type: none"> <li>• Observe the digital temperature/humidity display routinely and respond to humidity alarms.</li> <li>• Keep the sample chamber and accessories clean and dry.</li> <li>• Regenerate or replace the reusable stainless-steel desiccant before it is exhausted.</li> <li>• Store the instrument sealed if it will be idle for an extended period.</li> </ul> | <p><b>Routine care</b></p> <ul style="list-style-type: none"> <li>• Use trained personnel and follow the operating manual.</li> <li>• Clean KBr windows, ATR crystals and sample tools promptly after measurement.</li> <li>• Protect the source, detector, beam splitter and windows from prolonged moisture exposure.</li> <li>• Contact technical support if abnormal performance is observed.</li> </ul> |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

Specifications and configurations are subject to change without notice. Images are illustrative. Rev. A | 2026-08

## CONTACT / SUPPORT

### SpeSense Technologies Inc.

Add: No. 9 Yuhe Road, Guitang Subdistrict, Yuhua District, Changsha City

Tel: +86-15386478166

Email: [info@spesense.com](mailto:info@spesense.com)

Web: [www.spesense.com](http://www.spesense.com)



WhatsApp / WeChat